

DYM

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:76kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DYM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DMC; SMC;

Catalog #:A8222

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54808

Isotype:IgG

Swiss Prot:Q7RTS9

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein which is necessary for normal skeletal development and brain function. Mutations in this gene are associated with two types of recessive osteochondrodysplasia, Dyggve-Melchior-Clausen (DMC) dysplasia and Smith-McCort (SMC) dysplasia, which involve both skeletal defects and mental retardation.

To place an order, please [Click HERE](#).