

RNF31

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:120kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RNF31

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HOIP; ZIBRA;

Background:

The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions.

To place an order, please [Click HERE](#).

Catalog #:A8227

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55072

Isotype:IgG

Swiss Prot:Q96EP0

Purity:Affinity purification

For research use only.