

KANSL3

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:96kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KANSL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NSL3; Rcd1; KIAA1310;

To place an order, please [Click HERE](#).

Catalog #:A8234

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55683

Isotype:IgG

Swiss Prot:Q9P2N6

Purity:Affinity purification

For research use only.