

HPSE2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:67kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HPSE2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UFS; HPA2; HPR2; UFS1;

Catalog #:A8248

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:60495

Isotype:IgG

Swiss Prot:Q8WWQ2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a heparanase enzyme. The encoded protein is a endoglycosidase that degrades heparin sulfate proteoglycans located on the extracellular matrix and cell surface. This protein may be involved in biological processes involving remodeling of the extracellular matrix including angiogenesis and tumor progression. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).