

TMBIM1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:35kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human TMBIM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LFG3; RECS1; MST100; PP1201; MSTP100;

To place an order, please [Click HERE](#).

Catalog #:A8249

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:64114

Isotype:IgG

Swiss Prot:Q969X1

Purity:Affinity purification

For research use only.