

## DHX58

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**77kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human DHX58

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

LGP2; RLR-3; D11LGP2; D11lgp2e;

*To place an order, please [Click HERE](#).*

**Catalog #:**A8257

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**79132

**Isotype:**IgG

**Swiss Prot:**Q96C10

**Purity:**Affinity purification

For research use only.