

HAUS3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 70kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HAUS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IT1; dgt3; C4orf15;

Catalog #: A8258

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79441

Isotype: IgG

Swiss Prot: Q68CZ6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component of the HAUS augmin-like protein complex, which plays a key role in cytokinesis and mitosis. Disruption of the encoded protein causes mitotic defects resulting from fragmentation of centrosomes and microtubule destabilization. This gene shares its 5' exons with some transcripts from overlapping GeneID: 353497, which encodes a DNA polymerase. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).