

TCTN2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 77kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TCTN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MKS8; TECT2; C12orf38;

Catalog #: A8260

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79867

Isotype: IgG

Swiss Prot: Q96GX1

Purity: Affinity purification

For research use only.

Background:

This gene encodes a type I membrane protein that belongs to the tectonic family. Studies in mice suggest that this protein may be involved in hedgehog signaling, and essential for ciliogenesis.

Mutations in this gene are associated with Meckel syndrome type 8. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).