

ARL6

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ARL6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BBS3; RP55;

Catalog #: A8269

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84100

Isotype: IgG

Swiss Prot: Q9H0F7

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the ARF-like (ADP ribosylation factor-like) sub-family of the ARF family of GTP-binding proteins which are involved in regulation of intracellular traffic. Mutations in this gene are associated with Bardet-Biedl syndrome (BBS). Alternative splicing results in multiple transcript variants of this gene. A vision-specific transcript encoding a different protein has been described (PMID: 20333246).

To place an order, please [Click HERE](#).