

TP53INP1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human TP53INP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SIP; Teap; p53DINP1; TP53DINP1; TP53INP1A; TP53INP1B;

To place an order, please [Click HERE](#).

Catalog #:A8274

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:94241

Isotype:IgG

Swiss Prot:Q96A56

Purity:Affinity purification

For research use only.