

HYLS1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HYLS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HLS;

Catalog #:A8285

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:219844

Isotype:IgG

Swiss Prot:Q96M11

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein localized to the cytoplasm. Mutations in this gene are associated with hydroletharus syndrome. Multiple alternatively spliced variants, encoding the same protein, have been identified.

To place an order, please [Click HERE](#).