

## HYLS1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**34kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human HYLS1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HLS;

**Background:**

This gene encodes a protein localized to the cytoplasm. Mutations in this gene are associated with hydrolethalmus syndrome. Multiple alternatively spliced variants, encoding the same protein, have been identified.

*To place an order, please [Click HERE](#).*

**Catalog #:**A8285

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**219844

**Isotype:**IgG

**Swiss Prot:**Q96M11

**Purity:**Affinity purification

For research use only.