

SIK2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:104kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SIK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

QIK; SIK-2; SNF1LK2; LOH11CR11;

To place an order, please [Click HERE](#).

Catalog #:A8321

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23235

Isotype:IgG

Swiss Prot:Q9H0K1

Purity:Affinity purification

For research use only.