

## DHTKD1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 103kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human DHTKD1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CMT2Q; AMOXAD;

**Catalog #:** A8369

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 55526

**Isotype:** IgG

**Swiss Prot:** Q96HY7

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a component of a mitochondrial 2-oxoglutarate-dehydrogenase-complex-like protein involved in the degradation pathways of several amino acids, including lysine. Mutations in this gene are associated with 2-aminoadipic 2-oxoadipic aciduria and Charcot-Marie-Tooth Disease Type 2Q.

*To place an order, please [Click HERE](#).*