

PALB2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:131kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PALB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FANCN; PNCA3;

Catalog #:A8373

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79728

Isotype:IgG

Swiss Prot:Q86YC2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein that may function in tumor suppression. This protein binds to and colocalizes with the breast cancer 2 early onset protein (BRCA2) in nuclear foci and likely permits the stable intranuclear localization and accumulation of BRCA2.

To place an order, please [Click HERE](#).