

## PALB2

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**131kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PALB2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FANCN; PNCA3;

**Catalog #:**A8373

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**79728

**Isotype:**IgG

**Swiss Prot:**Q86YC2

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein that may function in tumor suppression. This protein binds to and colocalizes with the breast cancer 2 early onset protein (BRCA2) in nuclear foci and likely permits the stable intranuclear localization and accumulation of BRCA2.

*To place an order, please [Click HERE](#).*