

## BTLA

---

**Reactivity:** Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:1000 - 1:4000

**Calculated MW:** 33kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human BTLA

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BTLA1; CD272;

**Catalog #:** A8377

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 151888

**Isotype:** IgG

**Swiss Prot:** Q7Z6A9

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.

*To place an order, please [Click HERE](#).*