

ATAD2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:159kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human ATAD2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ANCCA; ATAD2;DKFZp667N1320;PRO2000;L16

Catalog #:A8382

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:29028

Isotype:IgG

Swiss Prot:Q6PL18

Purity:Affinity purification

For research use only.

Background:

ATAD2 is one of them encoding a conserved factor harbouring an AAA type ATPase domain and a bromodomain. ATAD2 controls chromatin dynamics, genome transcriptional activities and apoptotic cell response. This protein has 2 isoforms produced by alternative splicing.

To place an order, please [Click HERE](#).