

RSPO3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:31kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human RSPO3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PWTSR; THSD2; CRISTIN1;

Background:

This gene belongs to the R-spondin family. The encoded protein plays a role in the regulation of Wnt (wingless-type MMTV integration site family)/beta-catenin and Wnt/planar cell polarity (PCP) signaling pathways, which are involved in development, cell growth and disease pathogenesis. Genome-wide association studies suggest a correlation of this gene with bone mineral density and risk of fracture. This gene may be involved in tumor development.

To place an order, please [Click HERE](#).

Catalog #:A8389

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84870

Isotype:IgG

Swiss Prot:Q9BXY4

Purity:Affinity purification

For research use only.