

## MEX3C

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**69kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human MEX3C

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

RKHD2; BM-013; MEX-3C; RNF194;

**Catalog #:**A8393

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**51320

**Isotype:**IgG

**Swiss Prot:**Q5U5Q3

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of a family of proteins with two K homology (KH) RNA-binding domains and a C-terminal RING-finger domain. The protein interacts with mRNA via the KH domains, and the protein shuttles between the nucleus and cytoplasm. Polymorphisms in this gene may contribute to hypertension.

*To place an order, please [Click HERE](#).*