

RAB24

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:23kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RAB24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A8394

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:53917

Isotype:IgG

Swiss Prot:Q969Q5

Purity:Affinity purification

For research use only.

Background:

RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Olkkonen et al., 1993)

To place an order, please [Click HERE](#).