

## ESRRB

---

**Reactivity:**Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**56kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ESRRB

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ERR2; ERRb; ESRL2; NR3B2; DFNB35;

**Catalog #:**A8416

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2103

**Isotype:**IgG

**Swiss Prot:**O95718

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development.

*To place an order, please [Click HERE](#).*