

## GTF3A

---

**Reactivity:** Human Mouse

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 42kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GTF3A

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

AP2; TFIIIA;

**Catalog #:** A8426

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2971

**Isotype:** IgG

**Swiss Prot:** Q92664

**Purity:** Affinity purification

For research use only.

**Background:**

The product of this gene is a zinc finger protein with nine Cys

*To place an order, please [Click HERE](#).*