

## TGFB11

---

**Reactivity:**Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**50kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human TGFB11

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HIC5; ARA55; HIC-5; TSC-5;

**Background:**

This gene encodes a coactivator of the androgen receptor, a transcription factor which is activated by androgen and has a key role in male sexual differentiation. The encoded protein is thought to regulate androgen receptor activity and may have a role to play in the treatment of prostate cancer. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*

**Catalog #:**A8459

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7041

**Isotype:**IgG

**Swiss Prot:**O43294

**Purity:**Affinity purification

For research use only.