

## TGFB3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:**47kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human TGFB3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ARVD; RNHF; ARVD1; TGF-beta3;

**Background:**

This gene encodes a member of the TGF-beta family of proteins. The encoded protein is secreted and is involved in embryogenesis and cell differentiation. Defects in this gene are a cause of familial arrhythmogenic right ventricular dysplasia 1.

*To place an order, please [Click HERE](#).*

**Catalog #:**A8460

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7043

**Isotype:**IgG

**Swiss Prot:**P10600

**Purity:**Affinity purification

For research use only.