

RPP30

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RPP30

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TSG15;

To place an order, please [Click HERE](#).

Catalog #:A8479

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10556

Isotype:IgG

Swiss Prot:P78346

Purity:Affinity purification

For research use only.