

PSIP1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PSIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

p52; p75; PAIP; DFS70; LEDGF; PSIP2;

To place an order, please [Click HERE](#).

Catalog #: A8483

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11168

Isotype: IgG

Swiss Prot: O75475

Purity: Affinity purification

For research use only.