

## SYNPO

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 99kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human SYNPO

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A8484

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 11346

**Isotype:** IgG

**Swiss Prot:** Q8N3V7

**Purity:** Affinity purification

For research use only.

**Background:**

Synaptopodin is an actin-associated protein that may play a role in actin-based cell shape and motility. The name synaptopodin derives from the protein's associations with postsynaptic densities and dendritic spines and with renal podocytes (Mundel et al., 1997)

**To place an order, please [Click HERE](#).**