

EMCN

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EMCN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EMCN2; MUC14;

Background:

EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix (Kinoshita et al., 2001)

To place an order, please [Click HERE](#).

Catalog #:A8494

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51705

Isotype:IgG

Swiss Prot:Q9ULC0

Purity:Affinity purification

For research use only.