

RIPK4

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 92kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RIPK4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DIK; PKK; PPS2; RIP4; ANKK2; NKRD3; ANKRD3;

Catalog #: A8495

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54101

Isotype: IgG

Swiss Prot: P57078

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a serine/threonine protein kinase that interacts with protein kinase C-delta. The encoded protein can also activate NFkappaB and is required for keratinocyte differentiation. This kinase undergoes autophosphorylation.

To place an order, please [Click HERE](#).