

SMEK1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:95kDa

Observed MW:Refer to figures

Immunogen:

Recombinant Protein of human SMEK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

smk1; FLFL1; PP4R3; SMEK1; smk-1; PP4R3A; MSTP033; KIAA2010;

To place an order, please [Click HERE](#).

Catalog #:A8500

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55671

Isotype:IgG

Swiss Prot:Q6IN85

Purity:Affinity purification

For research use only.