

EHMT1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 141kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EHMT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GLP; GLP1; KMT1D; FP13812; EUHMTASE1; Eu-HMTase1; bA188C12.1;

Catalog #: A8513

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79813

Isotype: IgG

Swiss Prot: Q9H9B1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a histone methyltransferase that is part of the E2F6 complex, which represses transcription. The encoded protein methylates the Lys-9 position of histone H3, which tags it for transcriptional repression. This protein may be involved in the silencing of MYC- and E2F-responsive genes and therefore could play a role in the G0/G1 cell cycle transition. Defects in this gene are a cause of chromosome 9q subtelomeric deletion syndrome (9q-syndrome, also known as Kleefstra syndrome). Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).