

C19orf40

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant Protein of human C19orf40

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C19orf40;

Catalog #: A8520

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 91442

Isotype: IgG

Swiss Prot: Q9BTP7

Purity: Affinity purification

For research use only.

Background:

FAAP24 is a component of the Fanconi anemia (FA) core complex (see MIM 227650), which plays a crucial role in DNA damage response (Ciccia et al., 2007)

To place an order, please [Click HERE](#).