

RMI2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 16kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RMI2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BLAP18; C16orf75;

Catalog #: A8523

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 116028

Isotype: IgG

Swiss Prot: Q96E14

Purity: Affinity purification

For research use only.

Background:

RMI2 is a component of the BLM (RECQL3; MIM 604610) complex, which plays a role in homologous recombination-dependent DNA repair and is essential for genome stability (Xu et al., 2008)

To place an order, please [Click HERE](#).