

NOX1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 65kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NOX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MOX1; NOH1; NOH-1; GP91-2

Catalog #: A8527

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27035

Isotype: IgG

Swiss Prot: Q9Y5S8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the NADPH oxidase family of enzymes responsible for the catalytic one-electron transfer of oxygen to generate superoxide or hydrogen peroxide.

Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).