

GLP1R

Reactivity:Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GLP1R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A8547

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2740

Isotype:IgG

Swiss Prot:P43220

Purity:Affinity purification

For research use only.

Background:

GLP1R (glucagon-like peptide 1 receptor) is a protein-coding gene. Diseases associated with GLP1R include fasting hypoglycemia, and type 2 diabetes mellitus. GO annotations related to this gene include transmembrane signaling receptor activity and glucagon receptor activity. An important paralog of this gene is GHRHR.

To place an order, please [Click HERE](#).