

SLC24A4

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SLC24A4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9364

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 123041

Isotype: IgG

Swiss Prot: Q8NFF2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the potassium-dependent sodium/calcium exchanger protein family. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jul 2010]

To place an order, please [Click HERE](#).