

P2RX5

Reactivity: Human

Tested applications: WB FC

Recommended Dilution: WB 1:200 - 1:2000 FC 1:20 - 1:50

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human P2RX5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9366

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5026

Isotype: IgG

Swiss Prot: Q93086

Purity: Affinity purification

For research use only.

Background:

The product of this gene belongs to the family of purinoceptors for ATP. This receptor functions as a ligand-gated ion channel. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream gene, TAX1BP3 (Tax1 binding protein 3). [provided by RefSeq, Mar 2011]

To place an order, please [Click HERE](#).