

GOLPH3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GOLPH3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9371

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64083

Isotype: IgG

Swiss Prot: Q9H4A6

Purity: Affinity purification

For research use only.

Background:

The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is a peripheral membrane protein of the Golgi stack and may have a regulatory role in Golgi trafficking. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of these variants has not been determined. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).