

STAG2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100

Calculated MW: 141kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human STAG2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9400

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10735

Isotype: IgG

Swiss Prot: Q8N3U4

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a subunit of the cohesin complex, which regulates the separation of sister chromatids during cell division. Targeted inactivation of this gene results in chromatid cohesion defects and aneuploidy, suggesting that genetic disruption of cohesin is a cause of aneuploidy in human cancer. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]

To place an order, please [Click HERE](#).