

POLR1A

Reactivity:Human

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:195kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human POLR1A

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A9432

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:25885

Isotype:IgG

Swiss Prot:O95602

Purity:Affinity purification

For research use only.