

SLC47A1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 62kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human SLC47A1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9444

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55244

Isotype: IgG

Swiss Prot: Q96FL8

Purity: Affinity purification

For research use only.

Background:

This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).