

## FZD5

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**65kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human FZD5

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A9469

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7855

**Isotype:**IgG

**Swiss Prot:**Q13467

**Purity:**Affinity purification

For research use only.

**Background:**

Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD5 protein is believed to be the receptor for the Wnt5A ligand. [provided by RefSeq, Jul 2008]

*To place an order, please [Click HERE](#).*