

DSG2

Reactivity:Human Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:122kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DSG2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A9480

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1829

Isotype:IgG

Swiss Prot:Q14126

Purity:Affinity purification

For research use only.

Background:

Desmosomes are cell-cell junctions between epithelial, myocardial, and certain other cell types. This gene product is a calcium-binding transmembrane glycoprotein component of desmosomes in vertebrate epithelial cells. Currently, three desmoglein subfamily members have been identified and all are members of the cadherin cell adhesion molecule superfamily. These desmoglein gene family members are located in a cluster on chromosome 18. This second family member is expressed in colon, colon carcinoma, and other simple and stratified epithelial-derived cell lines. Mutations in this gene have been associated with arrhythmogenic right ventricular dysplasia, familial, 10. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).