

FTCD

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FTCD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9484

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10841

Isotype: IgG

Swiss Prot: O95954

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a bifunctional enzyme that channels 1-carbon units from formiminoglutamate, a metabolite of the histidine degradation pathway, to the folate pool.

Mutations in this gene are associated with glutamate formiminotransferase deficiency.

Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Dec 2009]

To place an order, please [Click HERE](#).