

AHI1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:137kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human AHI1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A9493

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54806

Isotype:IgG

Swiss Prot:Q8N157

Purity:Affinity purification

For research use only.

Background:

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]

To place an order, please [Click HERE](#).