

RQCD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:50

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human RQCD1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9522

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9125

Isotype: IgG

Swiss Prot: Q92600

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the highly conserved RCD1 protein family. The encoded protein is a transcriptional cofactor and a core protein of the CCR4-NOT complex. It may be involved in signal transduction as well as retinoic acid-regulated cell differentiation and development.

Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2012]

To place an order, please [Click HERE](#).