

UQCC2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:15kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human UQCC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A9986

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84300

Isotype:IgG

Swiss Prot:Q9BRT2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a nucleoid protein localized to the mitochondria inner membrane. The encoded protein affects regulation of insulin secretion, mitochondrial ATP production, and myogenesis through modulation of mitochondrial respiratory chain activity. [provided by RefSeq, Oct 2012]

To place an order, please [Click HERE](#).