

PELI3

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:51kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PELI3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A9999

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:246330

Isotype:IgG

Swiss Prot:Q8N2H9

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a scaffold protein and an intermediate signaling protein in the innate immune response pathway. The encoded protein helps transmit the immune response signal from Toll-like receptors to IRAK1/TRAF6 complexes. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011]

To place an order, please [Click HERE](#).