

## -Tubulin

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IF

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant Protein

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**To place an order, please [Click HERE](#).**

**Catalog #:**AC021

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**203068

**Swiss Prot:**P07437

For research use only.